

KDM4B / JMJD2B Rabbit mAb


ABZ2732

Product Overview


Name	KDM4B / JMJD2B Rabbit mAb
Accession(Primary)	O94953
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 122 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:1000; FC-1:200-1:1000; IP-1:200-1:1000

Target information (O94953)


Gene Name	KDM4B
About Gene Symbol	



Western blot analysis of KDM4B in 293, COS7 lysates using KDM4B antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using KDM4B/JMJD2B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of
KDM4B/JMJD2B(green) in Hela using KDM4B/JMJD2B
antibody,and DAPI(blue)