

ATG12 Rabbit mAb

ABZ2735


Product Overview

Name	ATG12 Rabbit mAb
Accession(Primary)	O94817
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (O94817)

Gene Name ATG12

About Gene Symbol



Western blot analysis of Apg12 in K562, 3T3, HeLa lysates using ATG12 antibody.

