

TPPP Rabbit mAb


ABZ2736

Product Overview



Name	TPPP Rabbit mAb
Accession(Primary)	O94811
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:100-1:200; FC-1:20-1:50

Target information (O94811)


Gene Name	TPPP
About Gene Symbol	




Immunocytochemistry analysis of TPPP (green) in LNCap using TPPP antibody, and DAPI (blue).



Western blot analysis of TPPP in mouse brain lysates using TPPP antibody.



Western blot analysis of TPPP in rat Brain lysates using TPPP antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Tubulin Polymerization Promoting Protein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.