

ALDH1L1 Rabbit mAb

ABZ2739


Product Overview

Name	ALDH1L1 Rabbit mAb
Accession(Primary)	O75891
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 99 kDa; Observed MW: 99 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (O75891)

Gene Name ALDH1L1

About Gene Symbol



Western blot analysis of ALDH1L1 in rat Brain, 3T3 lysates using ALDH1L1 antibody.

