

IDH1 Rabbit mAb


ABZ2740

Product Overview


Name	IDH1 Rabbit mAb
Accession(Primary)	O75874
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:10-1:100

Target information (O75874)

Gene Name	IDH1
About Gene Symbol	IDH1 (Isocitrate Dehydrogenase (NADP(+)) 1) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms



Western blot analysis of IDH1 in rat Brain, C6, 3T3, Hela lysates using IDH1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Isocitrate dehydrogenase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.