

BANF1/BAF Rabbit mAb


ABZ2741

Product Overview

Name	BANF1/BAF Rabbit mAb
Accession(Primary)	O75531
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 10 kDa; Observed MW: 10 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:100-1:200; IF-1:50-1:100; FC-1:50-1:100

Target information (O75531)

Gene Name	BANF1
About Gene Symbol	



Western blot analysis of BANF1 in Jurkat lysates
using BANF1 antibody.