

Fas Rabbit mAb

ABZ2742

Product Overview


Name	Fas Rabbit mAb
Accession(Primary)	P25445
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 40-50 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500

Target information (P25445)


Gene Name	FAS
About Gene Symbol	



Immunocytochemistry analysis of Fas (green) in HL-60 using Fas antibody, and DAPI (blue).



Western blot analysis of Fas in HeLa, Jurkat lysates using Fas antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using Fas antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.