

KIAA0652/ATG13 Rabbit mAb

ABZ2745


Product Overview

Name	KIAA0652/ATG13 Rabbit mAb
Accession(Primary)	O75143
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:50-1:100; IP-1:20-1:50


Target information (O75143)

Gene Name ATG13

About Gene Symbol



Western blot analysis of ATG13 in Raw264.7
lysates using ATG13 antibody.



Western blot analysis of ATG13 in K562, C6, HeLa lysates using ATG13 antibody.