

Syndecan 3 Rabbit mAb

ABZ2747


Product Overview

Name	Syndecan 3 Rabbit mAb
Accession(Primary)	O75056
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O75056)

Gene Name SDC3

About Gene Symbol



Western blot analysis of Syndecan 3 in K562, C6, 3T3 lysates using Syndecan 3 antibody.

