

JAK2 Rabbit mAb


ABZ2748

Product Overview

Name	JAK2 Rabbit mAb
Accession(Primary)	O60674
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 131 kDa; Observed MW: 131 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; IP-1:50-1:100

Target information (O60674)

Gene Name	JAK2
About Gene Symbol	



Western blot analysis of JAK2 in K562, 3T3 lysates using JAK2 antibody.