

KDM1/LSD1 Rabbit mAb


ABZ2749

Product Overview


Name	KDM1/LSD1 Rabbit mAb
Accession(Primary)	O60341
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 93 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:50-1:200; ICC/IF-1:100-1:200; FC-1:20; IP-1:20-1:50

Target information (O60341)


Gene Name	KDM1A
About Gene Symbol	



Western blot analysis of KDM1/LSD1 in Jurkat lysates using KDM1A antibody.



Western blot analysis of KDM1/LSD1 in HeLa, CHO-K1, C6, Jurkat lysates using KDM1/LSD1 antibody.



Immunocytochemistry analysis of
KDM1/LSD1(green) in HeLa using KDM1/LSD1
antibody, and DAPI(blue)