

TAK1 Rabbit mAb


ABZ2751

Product Overview

Name	TAK1 Rabbit mAb
Accession(Primary)	O43318
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:20-1:50; ICC/IF-1:200-1:1000; FC-1:20-1:50

Target information (O43318)

Gene Name	MAP3K7
About Gene Symbol	MAP3K7 (Mitogen-Activated Protein Kinase Kinase Kinase 7) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and TGF-β signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TAK1 in HeLa, A431
lysates using TAK1 antibody.