

Thioredoxin / TRX Rabbit mAb


ABZ2753

Product Overview


Name	Thioredoxin / TRX Rabbit mAb
Accession(Primary)	P10599
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; IP-1:10-1:100; ICC/IF-1:50-1:100

Target information (P10599)


Gene Name	TXN
About Gene Symbol	



Immunocytochemistry analysis of Thioredoxin (green) in HeLa using Thioredoxin antibody, and DAPI (blue).



Western blot analysis of Thioredoxin in Jurkat, HeLa lysates using Thioredoxin antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Thioredoxin/TRX antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.