

Thioredoxin / TRX Rabbit mAb

ABZ2754


Product Overview


Name	Thioredoxin / TRX Rabbit mAb
Accession(Primary)	P10599
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P10599)

Gene Name TXN

About Gene Symbol





Western blot analysis of Thioredoxin/TRX in HeLa, A549, HL-60, U2OS, U87-MG lysates using Thioredoxin/TRX antibody.