

JNK1 Rabbit mAb

ABZ2755


Product Overview

Name	JNK1 Rabbit mAb
Accession(Primary)	Q91Y86
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 46,54 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q91Y86)

Gene Name Mapk8

About Gene Symbol



Western blot analysis of JNK1 in Jurkat, rat Brain, C6 lysates using JNK1 antibody.

