

ELOVL5 Rabbit mAb

ABZ2756


Product Overview

Name	ELOVL5 Rabbit mAb
Accession(Primary)	Q9NYP7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:100; FC-1:20-1:50


Target information (Q9NYP7)

Gene Name ELOVL5

About Gene Symbol



Western blot analysis of ELOVL5 in Jurkat, HeLa lysates using ELOVL5 antibody.



Immunocytochemistry analysis of ELOVL5(green)
in HeLa using ELOVL5 antibody, and DAPI(blue)