

## Doublecortin Rabbit mAb


ABZ2759

### Product Overview

|                             |                                                           |
|-----------------------------|-----------------------------------------------------------|
| <b>Name</b>                 | Doublecortin Rabbit mAb                                   |
| <b>Accession(Primary)</b>   | O43602                                                    |
| <b>Host Species</b>         | Rabbit                                                    |
| <b>Clonality</b>            | Monoclonal                                                |
| <b>Isotype</b>              | IgG                                                       |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                         |
| <b>Calculated MW</b>        | Calculated MW: 41 kDa; Observed MW: 41 kDa                |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:100-1:200;<br>ICC/IF-1:100-1:200 |

### Target information (O43602)

|                          |     |
|--------------------------|-----|
| <b>Gene Name</b>         | DCX |
| <b>About Gene Symbol</b> |     |



Western blot analysis of Doublecortin in rat Brain  
lysates using Doublecortin antibody.