

Matrin 3 Rabbit mAb

ABZ2760

Product Overview


Name	Matrin 3 Rabbit mAb
Accession(Primary)	P43243
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 95 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:200-1:1000; FC-1:100-1:200

Target information (P43243)


Gene Name	MATR3
About Gene Symbol	MATR3 (Matrin 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Matrin 3 (green) in HL-60 using Matrin 3 antibody, and DAPI(blue).



Western blot analysis of Matrin 3 in Jurkat, HeLa lysates using Matrin 3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using Matrin 3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.