


PD-L2 Rabbit mAb

ABZ2762

Product Overview

Name	PD-L2 Rabbit mAb
Accession(Primary)	Q9BQ51
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9BQ51)

Gene Name PDCD1LG2**About Gene Symbol**

Western blot analysis of PDL2 in K562, C6 lysates using PD L2 antibody.

