

GATA1 Rabbit mAb


ABZ2763

Product Overview



Name	GATA1 Rabbit mAb
Accession(Primary)	P15976
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; ChIP-1:20; IHC-1:50-1:100; ICC/IF-1:50-1:100

Target information (P15976)

Gene Name	GATA1
About Gene Symbol	



Immunocytochemistry analysis of GATA1 (green) in K562 using GATA1 antibody, and DAPI (blue).



Western blot analysis of GATA1 in K562, HeLa lysates using GATA1 antibody.