

TNFAIP8 Rabbit mAb


ABZ2764

Product Overview

Name	TNFAIP8 Rabbit mAb
Accession(Primary)	O95379
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 23 kDa; Observed MW: 21-23 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:50-1:100

Target information (O95379)

Gene Name	TNFAIP8
About Gene Symbol	TNFAIP8 (TNF Alpha Induced Protein 8) is a cytokine and immune signaling protein involved in cell migration and invasion. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TNFAIP8 in CEM, K562
lysates using TNFAIP8 antibody.