

Tryptophanyl tRNA synthetase Rabbit mAb


ABZ2765

Product Overview


Name	Tryptophanyl tRNA synthetase Rabbit mAb
Accession(Primary)	P23381
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:20-1:50; FC-1:200-1:500; IHC-1:50-1:100

Target information (P23381)


Gene Name	WARS1
About Gene Symbol	



Western blot analysis of Tryptophanyl tRNA Synthetase/WRS in Jurkat, HeLa lysates using Tryptophan tRNA Ligase antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Tryptophanyl tRNA Synthetase/WRS antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Tryptophanyl
tRNA Synthetase/WRS(green) in HeLa using
Tryptophanyl tRNA Synthetase/WRS antibody, and
DAPI(blue)