

IL-4 Rabbit mAb


ABZ2766

Product Overview

Name	IL-4 Rabbit mAb
Accession(Primary)	P05112
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (P05112)

Gene Name	IL4
About Gene Symbol	



Western blot analysis of IL-4 in GSTIL-4 (aa25153) lysates using IL-4 antibody.