

Thrombomodulin Rabbit mAb


ABZ2767

Product Overview


Name	Thrombomodulin Rabbit mAb
Accession(Primary)	P07204
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 60 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:1000; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:10-1:100

Target information (P07204)

Gene Name	THBD
About Gene Symbol	



Immunocytochemistry analysis of Thrombomodulin (green) in A549 using Thrombomodulin antibody, and DAPI (blue).



Western blot analysis of Thrombomodulin in Jurkat lysates using Thrombomodulin antibody.