

TAF3 Rabbit mAb


ABZ2768

Product Overview


Name	TAF3 Rabbit mAb
Accession(Primary)	Q5VWG9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:500

Target information (Q5VWG9)

Gene Name	TAF3
About Gene Symbol	



Western blot analysis of TAF3 in Jurkat lysates using TAF3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TAF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.