

ATP Synthase C Rabbit mAb

ABZ2770


Product Overview


Name	ATP Synthase C Rabbit mAb
Accession(Primary)	P05496
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100

Target information (P05496)

Gene Name ATP5MC1

About Gene Symbol





Immunocytochemistry analysis of ATP synthase C (green) in Hela using ATP synthase C antibody, and DAPI(blue)