

MGST1 Rabbit mAb


ABZ2771

Product Overview


Name	MGST1 Rabbit mAb
Accession(Primary)	P10620
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100; FC-1:10-1:100

Target information (P10620)

Gene Name	MGST1
About Gene Symbol	



Western blot analysis of Microsomal GST1 in rat Liver lysates using MGST1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Microsomal Glutathione Stransferase 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.