

ARFRP1/ARP Rabbit mAb

ABZ2772


Product Overview

Name	ARFRP1/ARP Rabbit mAb
Accession(Primary)	Q13795
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500

Target information (Q13795)

Gene Name ARFRP1

About Gene Symbol



Western blot analysis of ARFRP1/ARP in C6 lysates using ARFRP1 antibody.

