

Asparagine Synthetase Rabbit mAb

ABZ2774


Product Overview

Name	Asparagine Synthetase Rabbit mAb
Accession(Primary)	P08243
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (P08243)

Gene Name ASNS


About Gene Symbol



Western blot analysis of Asparagine Synthetase in rat pancreas lysates using Asparagine Synthetase antibody.



Western blot analysis of Asparagine Synthetase in K562, HeLa lysates using Asparagine Synthetase antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Asparagine synthetase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.