

ENDOGL1 Rabbit mAb


ABZ2778

Product Overview

Name	ENDOGL1 Rabbit mAb
Accession(Primary)	Q9Y2C4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q9Y2C4)

Gene Name	EXOG
About Gene Symbol	EXOG (Exo/Endonuclease G) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ENDOGL1/ENGL in K562
lysates using ENDOGL1 antibody.