

Hsp90 alpha Rabbit mAb

ABZ2779


Product Overview

Name	Hsp90 alpha Rabbit mAb
Accession(Primary)	P07900
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50

Target information (P07900)

Gene Name HSP90AA1

About Gene Symbol



Western blot analysis of HSP90 in K562, C6, HeLa lysates using Hsp90 alpha antibody.

