

Lipoamide Dehydrogenase Rabbit mAb


ABZ2780

Product Overview


Name	Lipoamide Dehydrogenase Rabbit mAb
Accession(Primary)	P09622
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (P09622)


Gene Name	DLD
About Gene Symbol	DLD (Dihydrolipoamide Dehydrogenase) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and id



Immunocytochemistry analysis of Lipoamide Dehydrogenase (green) in HeLa using Lipoamide Dehydrogenase antibody, and DAPI (blue).



Western blot analysis of Lipoamide Dehydrogenase in K562, rat Brain, C6, 3T3, HeLa lysates using Lipoamide Dehydrogenase antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using Lipoamide Dehydrogenase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.