

PAH Rabbit mAb


ABZ2781

Product Overview

Name	PAH Rabbit mAb
Accession(Primary)	P00439
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; FC-1:200-1:500

Target information (P00439)

Gene Name	PAH
About Gene Symbol	



Western blot analysis of PAH in 3T3, C6, rat Brain, K562 lysates using PAH antibody.