

PTF1A Rabbit mAb

ABZ2785


Product Overview

Name	PTF1A Rabbit mAb
Accession(Primary)	Q9QX98
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q9QX98)

Gene Name Ptf1A

About Gene Symbol



Western blot analysis of PTF1A in K562, C6 lysates using PTF1A antibody.

