

Neurofascin Rabbit mAb


ABZ2786

Product Overview


Name	Neurofascin Rabbit mAb
Accession(Primary)	Q810U3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 138 kDa; Observed MW: 186,155,75 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:1000-1:5000

Target information (Q810U3)

Gene Name	Nfasc
About Gene Symbol	



Western blot analysis of Neurofascin in Human placenta lysates using Neurofascin antibody.



Western blot analysis of Neurofascin in mouse Brain, rat Brain, COS7, HeLa lysates using Neurofascin antibody.