

PFKM Rabbit mAb


ABZ2787

Product Overview


Name	PFKM Rabbit mAb
Accession(Primary)	P08237
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:100-1:500

Target information (P08237)

Gene Name	PFKM
About Gene Symbol	PFKM (Phosphofructokinase, Muscle) is a protein phosphatase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Fructose 6 Phosphate Kinase in C6, 3T3, HeLa lysates using Fructose 6 Phosphate Kinase antibody.



Immunocytochemistry analysis of Fructose 6 Phosphate Kinase(green) in HeLa using Fructose 6 Phosphate Kinase antibody, and DAPI(blue)