

# Cytokeratin 4 Rabbit mAb


ABZ2789

## Product Overview

|                             |                                                                                          |
|-----------------------------|------------------------------------------------------------------------------------------|
| <b>Name</b>                 | Cytokeratin 4 Rabbit mAb                                                                 |
| <b>Accession(Primary)</b>   | P19013                                                                                   |
| <b>Host Species</b>         | Rabbit                                                                                   |
| <b>Clonality</b>            | Monoclonal                                                                               |
| <b>Isotype</b>              | IgG                                                                                      |
| <b>Species Reactivity</b>   | Human                                                                                    |
| <b>Calculated MW</b>        | Calculated MW: 56 kDa; Observed MW: 56 kDa                                               |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:50-1:100; FC-1:100-1:500; IP-1:10-1:100 |

## Target information (P19013)

|                          |                                                                                                                                                                                                                                                                                     |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | KRT4                                                                                                                                                                                                                                                                                |
| <b>About Gene Symbol</b> | KRT4 (Keratin 4) is a keratin cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Cytokeratin 4 in C6, HeLa lysates using Cytokeratin 4 antibody.