

Cytokeratin 16/K16 Rabbit mAb

ABZ2791


Product Overview

Name	Cytokeratin 16/K16 Rabbit mAb
Accession(Primary)	P08779
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500


Target information (P08779)

Gene Name KRT16

About Gene Symbol



Western blot analysis of Cytokeratin 16 in HeLa lysates using Cytokeratin 16 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Cytokeratin 16 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.