

PBR Rabbit mAb


ABZ2793

Product Overview


Name	PBR Rabbit mAb
Accession(Primary)	B1AH88
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 11 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:1000; ICC/IF-1:200-1:1000; FC-1:50-1:200; IP-1:10-1:100

Target information (B1AH88)


Gene Name	TSPO
About Gene Symbol	



Western blot analysis of PBR/TSPO in 3T3, HeLa lysates using PBR antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PBR antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of PBR/TSPO
(green) in Hela using PBR/TSPO antibody, and
DAPI(blue)