

GFAP Rabbit mAb


ABZ2794

Product Overview

Name	GFAP Rabbit mAb
Accession(Primary)	P14136
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC:1:50-1:100; ICC/IF-1:100-1:200; FC-1:100-1:600

Target information (P14136)

Gene Name	GFAP
About Gene Symbol	



Western blot analysis of GFAP in rat Brain, Hela lysates using GFAP antibody.