

Profilin 1 Rabbit mAb


ABZ2796

Product Overview


Name	Profilin 1 Rabbit mAb
Accession(Primary)	P07737
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100; FC-1:100-1:200

Target information (P07737)


Gene Name	PFN1
About Gene Symbol	



Immunocytochemistry analysis of Profilin 1 (green) in A549 using Profilin 1 antibody, and DAPI (blue).



Western blot analysis of Profilin 1 in C6, 3T3, HeLa lysates using Profilin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Profilin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.