


TGF beta 1 Rabbit mAb

ABZ2797


Product Overview

Name	TGF beta 1 Rabbit mAb
Accession(Primary)	P04202
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P04202)

Gene Name Tgfb1**About Gene Symbol**

Western blot analysis of TGF beta 1 in Raw264.7 lysates using TGF beta 1 antibody.



Western blot analysis of TGF beta 1 in Jurkat, C6
lysates using TGF beta 1 antibody.