

Heme Oxygenase 1 Rabbit mAb


ABZ2798

Product Overview


Name	Heme Oxygenase 1 Rabbit mAb
Accession(Primary)	P09601
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:20

Target information (P09601)

Gene Name	HMOX1
About Gene Symbol	



Immunocytochemistry analysis of Heme Oxygenase 1 (green) in U87-MG using Heme Oxygenase 1 antibody, and DAPI (blue).



Western blot analysis of Heme Oxygenase 1 in 3T3 lysates using Heme Oxygenase 1 antibody.