

ITGA2B/CD41 Rabbit mAb


ABZ2800

Product Overview

Name	ITGA2B/CD41 Rabbit mAb
Accession(Primary)	P08514
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 113 kDa; Observed MW: 137 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (P08514)

Gene Name	ITGA2B
About Gene Symbol	



Western blot analysis of CD41 in C6, HeLa lysates using Integrin alpha IIb antibody.