

Rb Rabbit mAb


ABZ2801

Product Overview

Name	Rb Rabbit mAb
Accession(Primary)	P13405
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 105 kDa; Observed MW: 105 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:200-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (P13405)

Gene Name	Rb1
About Gene Symbol	



Western blot analysis of Rb in C6, Jurkat, HeLa lysates using Rb antibody.