


CD11b Rabbit mAb

ABZ2803

Product Overview

Name	CD11b Rabbit mAb
Accession(Primary)	P05555
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 127 kDa; Observed MW: 170 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:200-1:500; IP-1:20-1:50

Target information (P05555)

Gene Name Itgam**About Gene Symbol**

Western blot analysis of CD11b in Raw264.7
lysates using CD11b antibody.

