

c-Jun Rabbit mAb


ABZ2806

Product Overview


Name	c-Jun Rabbit mAb
Accession(Primary)	P05412
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 48,43 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:50-1:200; ICC/IF-1:20-1:50; FC-1:20; IP-1:20

Target information (P05412)


Gene Name	JUN
About Gene Symbol	JUN (Jun Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of c-Jun (green) in 293T using c-Jun antibody, and DAPI (blue).



Western blot analysis of c-Jun in 293, 293T lysates using c-Jun antibody.



Western blot analysis of c-Jun in 3T3 lysates using
c-Jun antibody