

ATF2 Rabbit mAb


ABZ2807

Product Overview


Name	ATF2 Rabbit mAb
Accession(Primary)	P15336
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 55 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; IP-1:20-1:50; ICC/IF-1:50-1:100

Target information (P15336)

Gene Name	ATF2
About Gene Symbol	



Immunocytochemistry analysis of ATF2 (green) in A549 using ATF2 antibody, and DAPI (blue).



Western blot analysis of ATF2 in Jurkat lysates using ATF2 antibody.