

14-3-3 alpha/beta Rabbit mAb


ABZ2809

Product Overview


Name	14-3-3 alpha/beta Rabbit mAb
Accession(Primary)	P31946
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (P31946)

Gene Name	YWHAB
About Gene Symbol	YWHAB is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of 41701 in K562, rat Brain, C6, 3T3, HeLa lysates using 14-3-3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using 14-3-3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.