

15 PGDH Rabbit mAb


ABZ2810

Product Overview

Name	15 PGDH Rabbit mAb
Accession(Primary)	P15428
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (P15428)

Gene Name	HPGD
About Gene Symbol	



Western blot analysis of 15PGDH in SW480
lysates using 15 PGDH antibody.